

---

## Technical

---

# Patient Risk Surveillance

## Product Compatibility Matrix LAB01795

### Description

The Product Compatibility Matrix outlines the software / hardware compatible with Patient Risk Surveillance (PRS). This document is intended to be referenced in relation to the 1.0.400 release of PRS.

### Product Affected

This bulletin outlines the critical software / hardware that are compatible with Patient Risk Surveillance:

- Smart Device Connectivity (SDC)
- Voalte Family
- Navicare Nursecall & Status board

### Issue Details

The product compatibility matrices contained in this bulletin define the critical software/hardware components that are compatible with Patient Risk Surveillance 1.0.400. All products/versions not listed in the matrices will require validation and an update to this technical bulletin.

## Browser Compatibility

The Patient Risk Surveillance shall be compatible with the following browsers:

- Chrome Version 76 or higher
- Microsoft Edge Version 84

## 3<sup>rd</sup> Party Integrations

Patient Risk Surveillance can be integrated to work with various 3<sup>rd</sup> party systems. The following matrixes outlines the 3<sup>rd</sup> party systems that are compatible with Patient Risk Surveillance 1.0.400.

Product	Compatible Versions
Smart Device Connectivity	1.2.400 and above
Voalte Server	3.7.17 and above
Navicare (including Nurse call & Status board)	3.9.600 and above

**COPYRIGHT © HILL-ROM SERVICES, INC., 2021. All rights reserved. No part of this text shall be reproduced or transmitted in any form or by any means, electronic or mechanical, including photocopying, recording, or by any information or retrieval system without written permission from HILL-ROM SERVICES, INC. (Hillrom). VersaCare, TotalCare, NaviCare, COMLinx, SideCom and Hill-rom are trademarks or registered trade-marks of Hill-rom Services, Inc. Windows is a registered trademark of Microsoft Corporation. All other products or brand names are trademarks or registered trademarks of their respective companies.**

The information contained in this guide is subject to change without notice. The only product warranty intended by Hill-rom is the express, written warranty accompanying the bill of sale to the original purchaser. Hill-rom makes no other warranty, express or implied, and in particular, makes no warranty of merchantability or fitness for any particular purpose. Additional hard copies of this guide can be obtained from Hill-rom by addressing your request to: Hill-rom Technical Support Department Telephone: (800) 445-3720.